

Initial Performance of the RatCAP, a PET Camera for Conscious Rat Brain Imaging


P. Vaska, C. Woody, D. Schlyer, V. Radeka, P. O'Connor, J.-F. Pratte, S. Shokouhi, S. Stoll, S. Junnarkar, M. Purschke, S.-J. Park, S. Southekal, V. Dzhordzhadze, W. Schiffer, D. Marsteller, D. Lee, A. Villanueva, S. Boose, A. Kandasamy, B. Yu, A. Kriplani, S. Krishnamoorthy

Brookhaven National Laboratory

J. Neill, M. Murphy, T. Aubele, R. Kristiansen
Long Island University

R. Lecomte and R. Fontaine
University of Sherbrooke

February 8, 2006



Animal Neuroscience with PET

Eliminate anesthesia:


- Why?
 - 1) removes confounding effects on neurochemistry
 - 2) enables stimulation during PET scanning and correlation w/ behavior
- How?
 - 1) miniaturize scanner
 - 2) attach directly to skull



Rat Conscious Animal PET


■ Mechanical

- FOV: 38 mm dia. x 18 mm axial
- OD: 72 mm
- 194 g counterweighted in rat-turn bowl



■ Detectors

- 12 blocks of 4x8 crystals = 384 crystals
- LSO: 2.2 mm x 2.2 mm x 5 mm long
- APD: Hamamatsu S8550




■ Electronics

- rigid-flex circuit board w/ custom ASICs
- small cable for all interconnections
- downstream: TSPM, VME daq, PC

System Features

- Singles mode
- No ADCs
- Custom ASIC CMOS 0.18 μ m
 - 32 preamps, bipolar shapers, zero-cross discriminators
 - ENC <1000 e⁻ rms, total power ~1 W
 - 32:1 multiplexer to single data line
 - 100 MHz clock: <100 ns deadtime per event per block
- TDC & Signal Processing Module (TSPM)
 - TDC with 1.1 ns time bins
- Data acquisition
 - VME FIFO & pdaq software on linux box
 - > 10 MB/sec w/ no deadtime



System Setup

■ Assembly

- >98% of detectors working

■ Gain variations

- ~ 2:1 max:min within block
- HV trimmed to reduce variations among blocks >> not optimized!

System Setup

■ Assembly


- >98% of detectors working

■ Gain variations

- ~ 2:1 max:min within block
- HV trimmed to reduce variations among blocks >> not optimized!

■ HV & energy threshold

- energy spectrum via threshold scan
- HV set to use full dynamic range
- common LLD voltage



System Setup

■ Assembly

- >98% of detectors working


■ Gain variations

- ~ 2:1 max:min within block
- HV trimmed to reduce variations among blocks >> not optimized!

■ HV & energy threshold

- energy spectrum via threshold scan
- HV set to use full dynamic range
- common LLD voltage

Avg LLD =
146 keV



System Setup

■ Assembly

- >98% of detectors working

■ Gain variations

- ~ 2:1 max:min within block
- HV trimmed to reduce variations among blocks >> not optimized!


Avg LLD =
146 keV

■ HV & energy threshold

- energy spectrum via threshold scan
- HV set to use full dynamic range
- common LLD voltage

■ Time window (offline)

- spectrum for each crystal pair
- 13.9 ns FWHM
- shifts ~10 ns not yet corrected
- >> using $2\tau = 40$ ns



System Setup

■ Assembly

- >98% of detectors working

■ Gain variations

- ~ 2:1 max:min within block
- HV trimmed to reduce variations among blocks >> not optimized!


Avg LLD =
146 keV

■ HV & energy threshold

- energy spectrum via threshold scan
- HV set to use full dynamic range
- common LLD voltage


■ Time window (offline)

- spectrum for each crystal pair
- 13.9 ns FWHM
- shifts ~10 ns not yet corrected
- >> using $2\tau = 40$ ns




Data Processing


- Sinogram binning (offline)
 - fully uncompressed, 3D, span=1 sinograms
 - prompt & delayed
- Image reconstruction: MLEM
 - Monte Carlo precalculated system matrix p_{ij} (SimSET)
 - 1 mm³ voxels
 - decays per voxel >10x expected in animal study
 - ~20 sec/iteration
- Quantitative corrections
 - scatter, attenuation - fixed, from microCT
 - detector-pair efficiencies ε_i
 - “direct inversion” = ratio of measured to simulated sinograms of uniform phantom
 - randoms R_i from delayed sino



Animal Conditioning - Why?



Animal Conditioning



Mounting on Conditioned Animal




Wearing Full System



Preliminary Performance

- Energy resolution - scatter rejection
 - avg 23% FWHM




Preliminary Performance

- Energy resolution - scatter rejection

- avg 23% FWHM

- Time resolution - randoms rejection

- avg 14 ns FWHM



$$R = 2\tau s_1 s_2$$

Preliminary Performance


- Energy resolution - scatter rejection
 - avg 23% FWHM
- Time resolution - randoms rejection
 - avg 14 ns FWHM
- Point source sensitivity - trues only
 - 0.7% @ 150 keV (0.4% @ 400 keV)

small-animal PET sensitivities
(@ LLD = 250 keV)

orig microPET (Cherry)	0.56%
ATLAS	1.8%
microPET R4	4.4%
microPET P4	2.3%
microPET II (Cherry proto)	2.3%
microPET Focus 220	3.4%
microPET Focus 120	7.7%



Preliminary Performance

- Energy resolution - scatter rejection
 - avg 23% FWHM
- Time resolution - randoms rejection
 - avg 14 ns FWHM
- Point source sensitivity - trues only
 - 0.7% @ 150 keV (0.4% @ 400 keV)
- Deadtime - bias




Preliminary Performance

- Energy resolution - scatter rejection
 - avg 23% FWHM
- Time resolution - randoms rejection
 - avg 14 ns FWHM
- Point source sensitivity - trues only
 - 0.7% @ 150 keV (0.4% @ 400 keV)
- Deadtime - bias
- Natural background from ^{176}Lu decay
 - calculate 2600 dps for whole system
 - true coincidence rate = 80 cps



Noise Equivalent Countrate

- True rate alone not useful for estimating SNR
- Must account for added noise from background
- NEC = effective counting efficiency
- Counting statistics \gg Poisson




$$N = T + S + R$$

$$T = N - S - R$$

$$\sigma_T = \sqrt{\sigma_N^2 + \sigma_R^2} = \sqrt{N + kR} = \sqrt{T + S + (1 + k)R}$$

$$\frac{\sigma_{NEC}}{NEC} = \frac{1}{\sqrt{NEC}} \equiv \frac{\sigma_T}{T} = \frac{\sqrt{T + S + (1 + k)R}}{T}$$

$$NEC = \frac{T^2}{T + S + (1 + k)R}, \quad k = 0, 1$$



if $S \ll T$ and $S \ll R$,

$$NEC \approx \frac{T}{1 + (2)(R/T)}$$


Preliminary Performance

■ NEC Results

- phantom fills FOV, 19 cc
- peak NEC 10 kcps @ 4 $\mu\text{Ci}/\text{cc}$
- 20 ns window: NEC to 14 kcps
- microPET proto: 4 kcps @ 8 $\mu\text{Ci}/\text{cc}$
- microPET R4: 45 kcps @ 6 $\mu\text{Ci}/\text{cc}$



■ Limitations

- Difficult to estimate where real data falls
- Phantom only within FOV: no OOFOM background randoms!
- NEC ignores non-statistical aspects of image quality



Spatial Resolution



- ^{68}Ge point source
 - ~1 mm dia.
 - $r = 0\text{-}16 \text{ mm}$
- 2D FBP
 - sinogram arc correction by linear resampling
 - ramp filter
- 3D Monte Carlo MLEM
 - 50 iterations
- Note:
 - arc correction parameters to be optimized
 - point source size NOT deconvolved



First Phantom Data

■ Rat striatum phantom

- 3.4:1 ratio
- microPET R4
 - FBP (ramp)
- RatCAP
 - MLEM
 - 25 iterations
 - post-smoothing with 2 mm FWHM Gaussian



First Rat Image

■ Rat data

- 517 g, 802 μCi ^{18}F -FDG i.p.
- 45 min awake uptake, then chloral hydrate euthanasia

■ microPET R4 scan

- 10 min, LLD = 250 keV, $2\tau = 10$ ns
- 3D MLEM (MAP w/ $\beta=0$, 20 iterations)

■ RatCAP scan

- 33 minutes livetime over 150 min scan
 - equivalent to 1.9 X decays of R4 scan
- Monte Carlo-based 3D MLEM
 - 200 iterations
 - randoms correction
 - no efficiency correction (yet)



First Rat Image


microPET R4



RatCAP
(coreg using pmod)




overlay



RatCAP v. microPET


- Sensitivity lower but not easily improved
- Randoms fractions:
 - Realistic rates?
 - large rat, but FDG ~0.8 mCi >> high end of expected count rates for neuroreceptor studies
 - average activity ~1/2 of mPET activity (scanned later and for 3X longer)
 - cancels w/ 20 ns window
 - microPET data:
 - ~5 % randoms fraction within subject



$$NEC \approx \frac{T}{1 + (2)(R/T)}$$

Goals for Timing

- Effects of planned improvements:
 - individualized, higher LLDs
 - T ~ same, S lower, R lower due to lower singles rates
 - variance reduction in randoms should give factor of 2!
 - to 1st order, thicker crystals do NOT change fraction (2x length gives 2x singles = 4x randoms but also 4x trues)
- Goal
 - Reduce R/T at FDG brain rates to ~10%
 - Now, estimate ~ 50% average, or 25% w/ randoms variance reduction
 - >> ~2X better timing = 5 ns FWHM to allow 10 ns window




$$NEC \approx \frac{T}{1 + (2)(R/T)}$$

Developments Underway

- Hardware:
 - ASIC
 - Add individually adjustable thresholds and/or VGA
 - Add ESD pad protection
 - Future: improved timing & some energy info (time-over-threshold?)
 - Improve sensitivity w/ longer crystals: parallax/weight tradeoff
 - MRI-compatible version & testing
- Software
 - Correct for timing offsets & narrow coincidence window
 - Validate quantitative corrections, including scatter from outside the FOV
 - Reduce variance caused by randoms correction by calculating from singles rates (no coinc proc deadtime)
- Animal studies
 - Quantify stress levels w/ corticosterone & behavior measures
 - Full neuroreceptor study
 - F-18 fluoride bone scan

Thanks!



DOE OBER funding
29