


SPLAT-MS: Single Particle Laser Ablation Time-of-flight Mass Spectrometer



Alla Imre and Dan Imre

AKNOLEDGMENTS

Jian Wang	BNL	Klaus Muller	SUNYSB
Gunner Senum	BNL	Peter Imrich	SUNYSB
Jeong-Ho Han	BNL	Wei Zhu	SUNYSB
Logan Chieffo	SHU	Ray Mugno	SUNYSB
Tim Onasch	BNL		
Susan Oatis	SHU	Paul O'Connor	BNL
Mike Alexander	PNNL	Jack Fried	BNL

OUTLINE

- Single particle MS basics
 - SPLAT-MS performance
 - SPLAT-MS beyond size and composition
 - Single spore M-D characterization
-
- SpectraMiner data analysis and visualization
 - The interactive dendogram
 - Exploring particle data

Single Particle Properties Real-Time Measurements

Here today


- Size
- Composition
- Density (effective density)
- Water uptake

Coming soon


- Shape
- Optical – scattering, phase function, fluorescence
- Quantitative soot content

WHY?

Particles come internally and externally mixed




Externally mixed






Internally mixed



Single Particle Mass Spectrometry: Basics




INLET






DETECTION





Scattered light signal from PMT, generated by 80 nm particle




AERODYNAMIC SIZING



MASS SPECTRUM OF AMMONIUM SULFATE/SUCCINIC ACID PARTICLE



SPLAT-MS PHOTO



FIELD DEPLOYMENTS

TX 2000 Air Quality Study Houston TX
August-September 2000



Aerosol Characterization Experiment Cheju
Island Korea April 2001







PM2.5 Technology Assessment and
Characterization Study, NY July 2001



Mercedes A170 Engine Characterization,
NTRC Oak Ridge TN February 2003



4 Little Atmospheric Particles




ABLATION WILL NEVER BE QUANTITATIVE
&
WE WILL NEVER IDENTIFY ORGANICS BY NAME


WE CAN DO BETTER

SEPARATE IONIZATION AND EVAPORATION

SEPARATION OF EVAPORATION FROM IONIZATION






Ablation by UV laser only




CO₂ laser evaporation UV ionization

Application to Automobile Exhaust Emission




Size Distributions of Particles Observed by SPLAT-MS




Spectra of a Few Little Exhaust Particles

SOOT




Spectra of a Few Little Exhaust Particles

Lube Additives



Spectra of a Few Little Exhaust Particles


PAHs




- 30** NO
- 152** 178-C₂H₂
- 178** phenanthrene
- 192** C₁ Alkyl phenanthrene
- 202** Pyrene
- 216** Alkyl pyrene
- 234** C₄ Alkyl phenanthrene
- 252** Benzo pyrene

Size, Composition & Density

$$D_{\text{aero}}/D_{\text{DMA}} = \text{Density}$$




Density of 200nm Particles




Density of the Little Exhaust Particles

Lube Additives



Density of the Little Exhaust Particles


PAHs



30	NO
152	178-C ₂ H ₂
178	phenanthrene
192	C ₁ Alkyl phenanthrene
202	Pyrene
216	Alkyl pyrene
234	C ₄ Alkyl phenanthrene
252	Benzo pyrene


Density of the Little Exhaust Particles

SOOT




Density of the Little Exhaust Particles

SOOT




SOOT




Individual soot particles measure 10-60 nanometers wide; aggregates are 1,000 nm (1 micron) wide. (Photo from 1997's MSL-1 mission.)

Size Density Relationship for Soot Particles



Size Density Relationship for Soot Particles



BIOLOGICAL WARFARE

Detection of Biological Warfare agents
requires high **sensitivity** & **selectivity**

The ultimate detection limit

SINGLE SPORE

High selectivity through

M-D characterization

SENSITIVITY REQUIREMENT

Biological Agent	*Detection Goal
Botulinum Toxin	200 org./ cm ³ second
Yersinia Pestis	200 org./ cm ³ second
Coxiella Burnetii	200 org./ cm ³ second
Rift Valley Fever Virus	40 org./ cm ³ second
Bacillus Anthracis	40 org./ cm ³ second


*Based on a typical soldier breathing 1000 liters/hour

**ESTIMATED SENSITIVITY FOR
SINGLE PARTICLE MASS SPECTROMETER
~1 org./ cm³ second**



SINGLE SPORE ANALYSIS

- Maximum sensitivity
- Zero background
- Stepwise identification and discrimination on the basis of:
 - Size
 - Shape
 - Bio-fluorescence
 - Mass spectral signature
 - Antibody reactivity
 - Select purify and collect individual spores for microscopic, DNA, and other analysis


5-D CHARACTERIZATION



SHAPE & FLUORESCENCE




SCHEMATIC of a FIELD DEPLOYABLE INSTRUMENT




WHAT DO WE DO WITH ALL THE DATA?

- We can generate Gb of data each day.
- The data are detailed and if properly mined can provide insight into understanding aerosol sources and fate.
- So much information for such a little brain.
- We do not want to reduce (loose) the richness of the data.
- We cannot grow the brain
- We must develop tools that will make it possible to navigate through the data with ease.
- The tool must provide a multi-level view all the way from the individual particle to the entire data set with a road in-between.

The Dendrogram or Classification Tree



Composition Data with High Temporal Resolution



CONCLUSION

- ▶ Compact SPMS
- ▶ ~50 particles/sec
- ▶ 20nm to 3.5micron
- ▶ Size
- ▶ Composition by CI
- ▶ Optical properties
- ▶ Density
- ▶ Hygroscopicity
- ▶ Shape
- ▶ Fluorescence

- ▶ Real-time analysis and visualization
- ▶ Expert driven data classification
- ▶ Comprehensive data analysis and visualization of gases, MS, size distribution, etc.