


The limits of lithography and its relevance to nanotechnology (or: what will they build in our parking lot ?!)




John Warren

07/31/2002





- Clean rooms
- Vibration isolation of equipment
- Superior temperature, EMI isolation, and humidity control
- Common interaction areas
- Connected to NSLS and Instrumentation Division
- 78,500 sq. ft. lab and user space

What is Nanotechnology?



Substrate for Molecular Wires
(Stormer and Willet)



Ni precipitates on dislocation
in silicon (Warren)

Laboratory clusters

Electron Microscopy

High-resolution structural and chemical probes

Materials Synthesis

Bulk, thin film material synthesis capabilities

Nanopatterning

E-beam and Ion-beam writer, pattern transfer

Scanning Probe Microscopy

AFM, STM


Ultrafast Short Wavelength Source

Short wavelength-short pulse lasers


Nanocenter NSLS Beamlines

Small angle X-ray scattering and microprobe


Nanopatterning Laboratory: Current Projects




Nanotemplate Directed Assembly of Soft Matter and Biomaterials




Charge injection and Transport in Nanoscale Materials: C. Creutz et al., Chemistry Dept., BNL




Hybrid Semiconductor-Superconductor Nanostructures: E. Mendez and F. Camino, Dept. of Physics and Astronomy, SUNYSB



Fresnel Zone Plate for X-Ray Microscopy at NSLS: C. Jacobson, Dept. of Physics and Astronomy, SUNYSB





A



B


SEM micrographs of central gap in nanoelectrode microstructure. The nanoelectrode array (A) is composed of SU-8, a UV-sensitive negative resist. The electrodes are made conductive by using directional vacuum evaporation to coat the top surface of the electrode with a conducting Au/Cr layer. As shown in (B), metal is not deposited on the vertical sidewalls of the electrode, and electron isolation is maintained between the two halves of the electrode and the substrate.

How will nanotechnology affect Instrumentation?




Nanoembossing: Is optical lithography obsolete?


mold in contact w. substrate




substrate melts < 250 ns




excimer laser radiation



substrate solidifies $t > 250$ ns





mold separated fm. substrate






Ultrafast and direct imprint of nanostructures in silicon
S. Y. Chou, C. K. Kelmel & J. Gu, Nature, vol. 417, 20 June 2002

Etched quartz embossing disc and silicon wafer after nanoembossing process:



Ultrafast and direct imprint of nanostructures in silicon
S. Y. Chou, C. K. Kuo & J. Gu, *nature*, vol. 417, 20 June 2002

What is the distinction between direct-write and replication methods?



Direct-write Technologies for Rapid Prototyping Applications


Throughput:

Exposure time per "pixel" =

"Sensitivity" / Intensity


Dip-Pen Nanolithography

* as small as 15 nm linewidths and ~ 5 nm spatial resolution



D. Piner, J. Zhu, F. Xu and S. Hong, C. A. Mirkin, "Dip-Pen Nanolithography", Science, 1999, 283, 661-63.


Resolution vs. Areal Throughput (1 sq cm = 10^8 um^2) !



Nanopatterning Lab: Primary Instruments



Sequential Pattern Generation:

High Resolution Electron Beam Pattern Generator (JEOL 9300FS or Leica VB6-HR)



Focused Ion Beam Pattern Generation:


JEOL 9855S with 30 Kv Ga source and high resolution SEM capability for 200 mm wafers



Advantages & Disadvantages of Direct Write Nanofab Methods:



AFM Based	highest resolution, slowest thruput, unproven
Electron Beam	20 nm resolution, medium thruput, proven technology
Ion Beam	100 nm resolution, medium thruput, many substrates (not just resist), Ga ions only!
Laser	lower resolution but reduced cost: no vacuum, no 100 KeV power supplies, etc.

SEM image of patterned hole array in PMMA with silicon substrate patterned by 100 keV JEOL 9300 E-beam pattern generator





$$\lambda = h / \{2mVe(1 + eV/2mc^2)\}^{1/2}$$

If $\lambda = .037 \text{ \AA}$ for 100 keV electrons, why can we pattern (at best) 20 nm dots in PMMA? 20 nm is 5400 greater than λ ! Yet optical lens's are "diffraction limited" ?





Approx. size
of 0.1 nm dot-
the resolution
of a good TEM,
or the STEM in
Biology

Fraunhofer diffraction geometry for circular aperture



Definition of resolution and spot size



For small angles, $a/2 = D \alpha$

Solve for D and substitute:


Or: $y = \lambda / 2 \alpha$

So the focused spot diameter is:

$$d = 2y = \lambda / \alpha$$

For the first minimum:

$$y = \lambda D / a$$






Using the Raleigh resolution criterion:

$$d = 1.22 \lambda / a$$

Spot size as a function of semi-angle:

$$d = 1.22 \lambda / a$$





Flux in = Flux out!

(Langmuir theorem for any optical system), so:

$$b(\text{source}) \times \text{area} \times dW = b(\text{image}) \times \text{area} \times dW$$

By definition of a solid angle:



$$dW = p r^2 / R^2$$

But $r = Ra$ for small a

$$\text{So } dW = p a^2$$

$$\text{Flux} = b(\text{source}) \times \text{area} \times p a^2$$

in the image plane.



For an electron probe, the flux is replaced by the current I , and the flux expression: $b \times \text{area} \times dW$ becomes:


$$I = b \left(\pi d^2 / 4 \right) \left(\pi a^2 \right)$$

Now we solve for the spot size diameter:

$$d_g = \left(2 I^{1/2} / \pi b \right) \left(1 / a \right)$$

and we have defined the probe size using geometrical rather than physical optics

Electron beam probe size based on diffraction and geometrical optics considerations for a 100 keV electron gun




W thermal emission:
 $b = 1 \times 10^{-11}$ amps / sq.
 Angstrom



Field emission:




$b = 2 \times 10^{-8}$ amps / sq.
 Angstrom



$$d = 1.22 \lambda / a$$

$$d_g = (2 I^{1/2} / p b) (1 / a)$$


E field for focusing electrostatic lens and sample electron trajectories




Electromagnetic
lens

$$\mathbf{F} = q(\mathbf{E} + \mathbf{v} \times \mathbf{B})$$

Spherical aberration for electrostatic lens - external rays are brought to focus at shorter distance than paraxial rays



Chromatic aberration for electrostatic lens - less energetic electrons are brought to focus in shorter distance than ones with high energy



$$d_c = C_c a \, dE / E$$

Summing the contributions for diffraction, brightness, spherical aberration and chromatic aberration, we have:

$$d_i := 0.61 \cdot \frac{\lambda}{\alpha_i}$$


$$dgf_i := \left(\frac{2}{\pi} \cdot \sqrt{\frac{I}{\beta}} \right) \cdot \frac{1}{\alpha_i}$$




$$dc_i := Cc \cdot \alpha_i \cdot \frac{\delta E}{E}$$

$$ds_i := .30 \cdot Cs \cdot (\alpha_i)^3$$


$$D_i := \sqrt{(d_i)^2 + (dgf_i)^2 + (ds_i)^2 + (dc_i)^2}$$




The *area's* are summed to get total spot size:




Beam diameter as a function of semi-angle for high and low accelerating voltage




Electron beam diameter as a function of gun brightness




Probe diameter as a function of current for a fixed semi-angle:




Why is it so difficult to eliminate spherical aberration for both electrostatic and electromagnetic lens's?



Sub-nanometer spot size's are obtainable for field emission-based instruments, but the minimum feature size in resist is still much larger. WHY?




Electron-matter interactions:



Assumptions for Monte Carlo electron scattering:

- 1) Elastic (no energy loss) scattering (attraction between electron and nucleus; repulsion between electron and electron cloud resulting in angular path deviations) *completely determines the path taken by the electron.*
- 2) Inelastic scattering (energy loss) *takes place continuously along the path followed by the electron rather by discrete events* (inner shell ionization, etc)
- 3) Actual atom positions are ignored; matter is treated as a continuum.



The average distance traveled by the electron between elastic scattering events, is l , the mean free path:

$$l = A / N_a r_{SE}$$

l as a function of electron energy and atomic number:


Element	Z	100 keV	10 keV
Carbon	6	1310 Å	170 Å
Silicon	14	1112 Å	127 Å
Copper	29	297 Å	35 Å
Gold	79	89 Å	10 Å

In the Monte Carlo simulation, a random number generator determines the step length and the scattering angle for the electron path.


$$\cos f = 1 - \frac{2 a \text{RND}}{(1 + a - \text{RND})}$$

$$y = 2 p \text{ RND}$$

Probability of scattering angle exceeding a specified value. Most angles are small, but 50 % > 1.5°



The average step length is 1, but s is > 2.3 1 10% of the time.



Program sequence:*

1. Find the step length for a given material and energy.

$$s = -l \ln (\text{RND})$$

2. Find the angular deviation in terms of the direction cosines:

$$f = \cos^{-1} (2 a \text{ RND} / 1 + a - \text{RND})$$

$$y = 2 p \text{ RND}$$

where a is a function of the energy and atomic number.

3. Use a form of the “Bethe stopping power” eqn to compute energy lost:

$$dE/dS = -78500 * (Z/AE) \ln (1.166E / J)$$

where $S = rs$ and J is the mean ionization potential.

4. Using the direction cosines, compute the new position: xyz and assign the new energy to the electron.



*courtesy D. C. Joy, Oak Ridge Nat'l Lab

Scattering simulation in PMMA membranes



25 keV

100 keV

500 A




1000 A




Simulation of scattering of 25 keV and 100 keV electrons in bulk silicon substrate:

25 keV




100 keV




Simulation of 100 keV electrons with low Z and high Z elements


Silicon



Tungsten




Energy contours in 1000 Å PMMA film exposed by 100 keV electron beam



Summary:

- 1) The “limits” of resolution for any electron probe system are determined as much by electron-substrate interactions as the optics.
- 2) **Electron-resist and electron-substrate interactions must be optimized as we approach genuine nanoscale geometries.**

1



2



1. Sub-10 nm Electron Beam Lithography using Inorganic Resist, K. Yamamzaki et al. Proc. SPIE 3997 (2000) 458.
2. Univ. of Glasgow, unpublished rept.